

Method Procedure

Document N	lumber:	SASoM/METHOD/013.v4
Title:	Use of Fast protein norr	Green FCF for Reverse Phase Protein Array (RPPA) nalisation.
Version:	v4	
Author:	Peter Muller	1

Effective from:	21/03/2020	
Valid to:	20/03/2022	

SOP History		
Number	Date	Reason for Change
v1	03/04/2014	
v2	21/03/2016	Update
v3	21/03/2018	Update
v4	21/03/20	Biennial Update

1.0 Purpose –

This SOP describes the current procedure for using Fast Green FCF for Reverse Phase Protein Array (RPPA) protein normalisation in Laboratory 248 at the St Andrews School of Medicine (SASoM).

2.0 Scope -

This SOP applies to the staff in the SASoM involved in RPPA work in Laboratory 248 at the St Andrews School of Medicine (SASoM).

3.0 Responsibilities -

All staff must have read and signed the relevant risk assessment documents before performing this procedure.

4.0 Procedure –

Fast Green FCF dye is used for the general staining of proteins on porous nitrocellulose films for sample protein quantification. It is therefore particularly suited to Reverse Phase Protein Arrays (RPPA) where it can be used to quantify the 'total protein' deposited on nitrocellulose-coated 'Fastslides'.

Fast Green was purchased from Sigma (F7258-25G).



school of MEDICINE

Method Procedure

Be careful with the Fast Green though as it is hygroscopic so even a small amount left on the bench will attract water very quickly and the whole bench will look like St.Patrick's day in Dublin!

Make up 1 litre of Destaining Buffer as follows:

100ml Glacial Acetic Acid 300ml EtOH 600ml Distilled Water (green Elga tap)

Make up 50mL of Fast Green FCF dye (0.005% w/v) as follows:

2.5mg (0.0025g) of Fast Green FCF 50mL of Destaining Buffer Mix thoroughly until in suspension.

Procedure:

- 1. Rinse slides for 1 minute in PBS.
- 2. Incubate slides in 0.005% w/w Fast Green FCF dye solution for 2hrs at room temperature.
- 3. Rinse slides in Destaining buffer (3 x 5mins).
- 4. Incubate slides in Destaining buffer for 2hrs at room temperature, keeping well agitated throughout.
- 5. Air dry and read on the Licor Odyssey scanner. Images can be captured on both the 680nm (red) and 800nm (green) channels.

5.0 Personal protection -

A Howie laboratory coat and lab gloves must be worn at all times.

6.0 Spillages -

Always clean up any spills immediately after use, only you know what you have spill and are aware of its hazard.

Spillages should be mopped up with paper towel, disinfected with 70% ethanol and finally washed with distilled water and a damp cloth.

7.0 Related documents –

- 8.1 Risk assessments RA/GEN/009 (General Laboratory Safety)
- 8.1 SOPs SASoM/METHOD/038/039/040



Method Procedure

8.0 Approval and sign off –

Author:					
Name:	Peter Mullen				
Position:	Research Fellow				
Signature:		Date:			
Management Approval:					
Name:	Peter Mullen				
Position:	Research Fellow				
Signature:		Date:			
QA release by:					
Name:	Alex MacLellan				
Position:	QA Manager				
Signature:		Date:			





Method Procedure

STANDARD OPERATING PROCEDURE

Please sign below to indicate you have read this S.O.P and understand the procedures involved.

NAME	POSITION HELD	SIGNATURE	DATE