

Document Number:		SASoM/METHOD/037.v5
Title:	Preparation Spotting.	of 2-fold Sample Dilutions from FFPE for RPPA
Version:	V5	
Author:	Peter Muller	n

Effective from:	01/02/2021		
Valid to:	01/02/2023		

SOP History		
Number	Date	Reason for Change
V1	01/02/2013	Original
V2	01/02/2015	Update
V3	01/02/2017	Update
V4	01/02/2019	Update
V5	01/02/2021	Update

1.0 Purpose –

This SOP describes the current procedure for preparing serial dilutions of protein samples prior to spotting onto RPPA slides using the MGII Robotic Spotter in Laboratory 248 at the St Andrews School of Medicine (SASoM).

2.0 Scope -

This SOP applies to the staff in the SASoM involved in the preparation of 2-fold sample dilutions for RPPA spotting.

3.0 Responsibilities -

All staff involved with in the preparation of 2-fold sample dilutions from FFPE-derived protein samples for RPPA spotting are responsible for ensuring that the methods are followed in accordance with this SOP.

All staff must have read and signed the relevant risk assessment documents before performing this procedure.



4.0 Procedure –

Prepare protein lysates from FFPE blocks as described in the relevant SOP (SASoM-METHOD-069-Preparation of protein from FFPE (Paraffin) blocks).

Determine protein concentration of all samples using the Biorad RC-DC protein Assay kit (SASoM-METHOD-070-Protein Determination by the Biorad RC-DC Protein Assay). Keep samples on ice at all times.

A series of at least five 2-fold dilutions must be prepared for RPPA spotting onto Fast®slides in order that standard 4-parameter modelling with MicroVigene spot analysis software can be performed. (If 5-parameter modelling with MicroVigene is required, then at least six dilutions will be necessary).

For optimal spot recognition, sample dilutions should ideally start at 1mg/ml $1000\mu g/1000\mu l = 100\mu g/100\mu l$ including $25\mu l$ of RPPA denaturing buffer. All samples should therefore contain $100\mu g$ of total protein in a final volume of $100\mu l$, of which $25\mu l$ will be 4X Mercapthoethanol Reducing Buffer. The difference in volume from sample to sample is made up with fully supplemented Lysis Buffer (see above).

As a guideline, all sample lysates must therefore have a protein concentration greater than $133\mu g/100\mu l$ if 1mg/ml is to be adopted as the starting concentration for the serial dilutions.

Using small (0.6mL) centritubes, make up the first dilution of each sample with 400µg of total protein (differing volumes) and 50µl (1/4 volume) of 4X Mercapthoethanol Reducing Buffer. Make up the total volume of each sample to 200µl with fully-supplemented Lysis Buffer. Label each sample as 1A, 2A, 3A etc.

Denature all samples at 60°C for 60mins.

Make up serial dilutions as follows;

1A = $60 \mu l$ neat 1B = $120\mu l$ A + $120\mu l$ complete lysis buffer 1C = $100\mu l$ B + $100\mu l$ complete lysis buffer 1D = $80\mu l$ C + $80\mu l$ complete lysis buffer 1E = $60\mu l$ D + $60\mu l$ complete lysis buffer

Add 40µl of each sample dilution (1A-E, 2A-E, 3A-E etc) to appropriate wells of a 96well ROUND-BOTTOM tray, taking into consideration the positioning of the pins. ENSURE you mark which is to be the front and rear of the tray for when it gets loaded into the MGII Biobank.

Seal the tray with Parafilm and store at -80°C prior to spotting on the MGII Robotic spotter.

Trays can be placed back in the freezer after spotting to be subsequently re-spotted at a later date with no significant loss of activity.



5.0 Personal protection -

A Howie coat must be worn at all times. Gloves as specified in the appropriate COSHH RA.

6.0 Spillages –

Always clean up any spills immediately after use, only you know what you have spilt and are aware of its hazard.

Spillages should be mopped up with paper towel, disinfected with 70% ethanol and finally washed with detergent.

7.0 Training –

All staff should under go training in this technique before performing the procedure.

8.0 Related documents -

- 8.1 Risk assessments COSHH/028
- 8.2 SOP SASoM/METHOD/022 Protein Extraction from Human Tissues and Primary Xenografts
- 8.3 SASoM/METHOD/023 Protein Determination by the Bicinchoninic Acid (BCA) Assay.
- 8.4 SASoM/METHOD/006 Preparation of Protein Lysates from Adherent Cell Cultures



9.0 Approval and sign off –				
Author:				
Name:	Peter Mullen			
Position:	Research Fellow			
Signature:	Peter Muller	Date: 11/02/2021		
Management App	roval:			
Name:	Peter Mullen			
Position:	Research Fellow			
Signature:	Peter Muller	Date: 11/02/2021		
QA release by:				
Name:	Alex MacLellan			
Position:	QA Manager			
Signature:	Cales markeller	Date: 11/02/2021		

 \sim



STANDARD OPERATING PROCEDURE

Please sign below to indicate you have read this S.O.P and understand the procedures involved.

NAME	POSITION HELD	SIGNATURE	DATE
		•	