



Document Number: SASoM/METHOD/046.v5

Title: Cloning and validating a 19mer shRNA hairpin in pLKO.1-puro.

Version: v5

Author: Paul Reynolds

Effective from:	10/08/2021
Valid to:	09/08/2023

SOP History		
Number	Date	Reason for Change
v1	29/08/2013	Original
v2	10/08/2015	Update
v3	10/08/2017	Update
V4	10/08/2019	Update
V5	10/08/21	Update

1.0 Purpose –

This SOP describes the current procedure for cloning a 19mer shRNA hairpin to pLKO.1puro for use in Laboratory 248 at the St Andrews School of Medicine (SASoM).

2.0 Scope –

This SOP applies to the staff in the SASoM involved in cloning and validating a 19mer shRNA hairpin to pLKO.1-puro.

3.0 Responsibilities –

All staff involved in cell culture are responsible for ensuring that the methods are followed in accordance with this SOP.

All staff must have read and signed the relevant risk assessment documents before performing this procedure.



4.0 Procedure –

Steps:

Digest/CIP treat pLKO.1-puro empty vector
Design shRNA oligos
Anneal oligos
Ligate oligos and vector
Grow and purify DNA
Test digest and Sequence insert
Make virus, infect cells and perform Western blotting.

Digest/CIP Treat Vector: pLKO.1-puro

1 µg pLKO.1-puro vector maxiprep DNA
2 µL NEB buffer 1
1 µL AgeI (NEB)
1 µL EcoRI (NEB)
H₂O to 20 µL
Add 0.5 µL CIP
Incubate 37° 1 hr

Run digest plus loading dye in an agarose gel and extract DNA band using a Gel purification kit (Qiagen). Elute in 50 µL of MilliQ water, do NOT follow the elution instructions from the kit. (see protocols for agarose gels).

Design Short-hairpin Oligos:

Find 19-mer sequence to target gene

Tip—Use sequences that are verified in publications and have pictures of western blot knockdown at protein level.

e.g. Use the program psicoligomaker 1.5 from Tyler Jacks lab to design oligos

<http://web.mit.edu/ccr/labs/jacks/protocols/pSico.html>

Directions below:

For pLKO.1:

Type 19-mer gene-specific sequence in “Target”

Click “Settings”



Method Procedure

Change end sequences

loop = TTCAAGAGA

sense 5' = CCGG

antisense 3' =

sense 3' = G

antisense 5' = AATTC

Click "Done"

Click "Make oligos"

The "For Oligo" is your forward oligo

The "Rev Oligo" is your reverse oligo

Order the oligos from Invitrogen, using the following criteria:

Desalted, 50nmol, 5'phosphorylated.

Anneal Oligos:

Centrifuge tubes briefly, 10 sec at 10,000rpm. Resuspend oligos at 1nmol/μL with sterile MilliQ water.

Set up water bath at precisely 70°C. Check with thermometer.

Annealing solution: in an eppendorf tube,

1 μL of diluted forward oligo (1:10 dilution in MilliQ water of stock)

1 μL of diluted reverse oligo (1:10 dilution in MilliQ water of stock)

5 μL of 10x NEB buffer 2

43 μL MilliQ water to a final volume of 50μL.

Incubate 4 min at 95°C, then

Incubate 10 min at precisely 70°C in water bath - use a thermometer (in water bath or 1L beaker of H₂O heated by a Bunsen) and allow the water to cool to RT (takes a few hours). Then move to 4°C.



Ligate Oligos and Vector:

Use Roche Rapid Ligation Kit:

- 5 μ L of water
- 1 μ L of annealed oligos (ctl no oligo and dilutions 1:10, 1:100, 1:1000)
- 2 μ L of digested, CIP treated, gel purified vector
- 2 μ L 5x DNA dilution buffer
- 10 μ L 2x ligation buffer
- 1 μ L ligase

Incubate at RT for 1-24 hrs

Grow and Purify DNA:

Transform DH5 α with ligation mix:

- Thaw DH5 α on ice (subcloning grade from Invitrogen)
- Add ligation mix to aliquoted 50 μ L DH5 α
- Incubate on ice 20 min
- Heat shock at 42°C for 30sec in heat block
- Incubate on ice 10 min
- Add transformed bugs to 500 μ L LB (no selection)
- Shake at 37°C for 60 min
- Plate 200 μ L on LB /100 μ g/mL carbenicillin 10cm bacterial plates
- Incubate, inverted at 37°C overnight in bacterial incubator
- Pick colonies into culture tubes:
- Grow up in 3 mL of LB (100 μ g/mL carbenicillin (Sigma))
- Mini-prep DNA using kit.

KEEP ALL CULTURES IN FRIDGE (even empty tubes) AND ALL DNA UNTIL CLONES ARE VERIFIED BY DIGESTION AND SEQUENCING.

Test Digest and Sequence Insert:

Test digest

For pLKO.1

- 7 μ L MilliQ water and 1ug of mini-prep DNA
- 1 μ L NEB buffer 1
- 0.1 μ L BSA
- 1 μ L AgeI (NEB)
- 1 μ L SpeI (NEB)

Incubate at 37°C 1 hr

Test on 1% agarose gel (follow agarose gel protocols)



Expected result

For pLKO.1-puro:

Loss of 469 bp band if insert is correct (note: oligos should destroy Agel site, however if sense 5' oligo starts CCGGT Agel site is not destroyed and there will be a 469 bp to 528 bp shift if the insert is correct).

Sequencing primer:

pLKO.1-seq: 5' GACTATCATATGCTTACCGT

Use this primer with DNA and follow protocol for sequencing at Dundee sequencing facility for those that pass the digestion test.

If sequence correct, make a glycerol stock from contents of bacterial culture:

750 µL of bacterial culture

250 µL of 80% sterile glycerol

in a cryovial. Mix and store at -80°C in Reynolds Lab glycerol stock box.

AT THIS STAGE YOU MAY DISCARD OTHER, NON-VALIDATED CLONE MATERIALS.

Culture validated bacterial stock to make maxiprep DNA (follow maxiprep DNA protocol).

Use this DNA to make virus and infect recipient cells (follow lentivirus protocols).

Only when protein of interest is knocked down in cells, as assessed by western blotting, can the shRNA be considered as validated.

5.0 Personal protection -

A Howie laboratory coat and lab gloves must be worn at all times.

6.0 Spillages -

Always clean up any spills immediately after use, only you know what you have spilt and are aware of its hazard.

Spillages should be mopped up with paper towel, disinfected with 70% ethanol and finally washed with 100% ethanol or water.

7.0 Training -

All staff should be trained in sterile TC techniques before starting any TC work



8.0 Related documents –

8.1 Risk assessments –

RA/BIOL/04 - Culture of primary and established cell lines

RA/GM/001 - Working with Category 2 (Cat 2) viruses

8.2 SOPs -

SASoM/METHOD/047 – siRNA transfection of mammalian cells

SASoM/METHOD/048 – DNA Sequencing at the University of Dundee

SASoM/METHOD/050 Biorad Mini Western Blotting

SASoM/METHOD/052 Agarose Gels

9.0 Approval and sign off –

Author:

Name: Paul Reynolds

Position: Principal Investigator

Signature:

Date: 07/10/2021

Management Approval:

Name: Peter Mullen

Position: SOP Administrator

Signature:

Date: 28/09/2021

QA release by:

Name: John O' Connor

Position: QA Manager

Signature:

Date: 28/09/2021

