

Method Procedure

Document I	umber: SASoM/METHOD/068.v3	
Title:	Preparation of 4X Denaturation buffer for RPPA	
Version:	v3	
Author:	Peter Mullen	

Effective from:	12/08/2018	
Valid to:	11/08/2020	

SOP History		
Number	Date	Reason for Change
v1	12/08/2014	Original
V2	12/08/2016	Biennial Review
V3	12/08/2018	Update

1.0 Purpose -

This SOP describes the current procedure for preparing 4X Denaturation buffer for RPPA in Laboratory 248 at the St Andrews School of Medicine (SASoM).

2.0 Scope -

This SOP applies to the staff in the SASoM involved in RPPA work in Lab 248.

3.0 Responsibilities -

All staff involved in cell culture are responsible for ensuring that the methods are followed in accordance with this SOP.

All staff must have read and signed the relevant risk assessment documents before performing this procedure.



Method Procedure

4.0 Procedure –

This SOP details how to make up 4X denaturation buffer for RPPA serial dilutions.

The composition of the (4X) buffer is as follows: 250mM Tris, pH6.8 (final concentration when diluted will be 62.5mM) 8% SDS (final concentration when diluted will be 2%) 35% glycerol (final concentration when diluted will be 8.75%) 10% Mercaptoethanol (final concentration when diluted will be 2.5%)

1. Make up a stock 1M Tris buffer by adding 12.114g of Tris base to a beaker. Add approximately 50mL of distilled water (green Elga taps) and correct to pH6.8 by adding HCI. As a result of the HCL being added, the temperature (and pH) of the buffer may increase. Allow the buffer to cool down before re-adjusting if necessary. Make up to a final volume of 100mL.



5.0 Personal protection -

A Howie laboratory coat and lab gloves must be worn at all times.

6.0 Spillages -

Always clean up any spills immediately after use, only you know what you have spilt and are aware of its hazard.

Spillages should be mopped up with paper towel, disinfected with 70% ethanol and finally washed with warm water.



Method Procedure

school of MEDICINE

7.0 Training -

8.0 Related documents -

8.1 Risk assessments – RA/COSHH/028 (RPPA)

9.0 Approval and sign off -

Author:							
Name:	Peter Mullen						
Position:	Research Fellow						
Signature:	Da	ate:					
Management Approval:							
Name:	Mary Wilson						
Position:	Laboratory Manager						
Signature:	Da	ate:					
QA release by:							
Name:	Alex MacLellan						
Position:	QA Manager						
Signature:	Da	ate:					