



Document Number: SASoM/METHOD/085.v3

Title: MCM2 Immunohistochemistry using Dako Link48

Version: v3

Author: In Hwa Um

Effective from:	15/03/20
Valid to:	14/03/22

SOP History		
Number	Date	Reason for Change
v1	15/03/2016	Original
v2	15/03/2018	Update
v3	15/03/2020	Update

1.0 Purpose –

This SOP describes the current procedure for staining cytology samples and formalin fixed paraffin embedded tissue sections with MCM2(Cytosystem) antibody using Dako Link 48 in Laboratory 248 at the St Andrews School of Medicine (SASoM).

2.0 Scope –

This SOP applies to the staff in the SASoM involved with MCM2 immunohistochemistry using Dako Link48.

3.0 Responsibilities –

All staff involved in immunohistochemistry using DAKO Link48 are responsible for ensuring that the methods are followed in accordance with this SOP. All staff must have read and signed the relevant risk assessment documents before performing this procedure.



4.0 Procedure –

MCM2 immunohistochemistry in Dako Link48

1. Turn on Dako automation system UPS and computer.
2. Double-click the DakoLink icon on the computer desktop and log in.



3. Insert slide information in New slides tab

DakoLink

System Administration Help

New slides Workflow Instruments Completed

New slides

Patient case

Case number: 4567 Accession number: Hospital: Downtown Hospital

First name: John Middle name: L Last name: Doe Accession date: 7/10/2012 Pathologist: Dr. A

Slides

Protocols: Melanosome HMB45 Blocks: Type: Preparations: FFPE Drop zones: Control slide: Slide notes: Quantity: 1 Add slide

- 3.1 Add Patient case info
 - 3.2 Protocol-choose primary antibody '**MCM2cyt**', which includes '**USTAN IHC DAB**' visualization system (Always double check if 'MCM2 cyt' has '**USTAN IHC DAB** visualization system linked.)
 - 3.3 Block-skip
 - 3.4 Type-skip
 - 3.5 Preparation-choose either FFPE or Cytology
 - 3.6 Drop zones-choose the centre for cytology and appropriate zones for FFPE
 - 3.7 Slide notes-if necessary
 - 3.8 Quantity
 - 3.9 Print and apply the slide labels
 - 3.10 Click 'Case complete'
 - 3.11 Stick the label on the slides
4. Put the slide(s) in non-metallic rack. (No dewaxing and rehydration steps required)



Manual antigen retrieval

- 4.1 Make up EDTA-Tris 1 L solution:
- 10mM Tris Base = 1.21g
1mM EDTA = 0.37g
Tween = 0.5ml
dH2O = 1L




Method Procedure

Add the EDTA and Tris Base to a 1L flask and add the dH2O and mix on magnetic stirrer. Ensure that the pH is 9. Add the Tween.

- 4.2 Pour the antigen retrieval solution into the microwave pressure cooker. Screw on the lid and put into the microwave on high power to heat up for 13mins.
- 4.3 Perform antigen retrieval on slides by putting slide rack into pre-heated pressure cooker. Screw on the lid and place the red weight over the vent. Put back into the microwave and heat on high for 5mins – there should be an audible hissing sound towards the end of the procedure.
- 4.3 Using heat resistant gloves and eye protection carefully lift the pressure cooker from within the microwave and place in sink. Place eye protection on. With the heat resistant gloves release the pressure by removing the weight from the top of the cooker. Stand well back immediately once the weight is removed. Do not stand over the pressure cooker during this step.
- 4.4 Allow the slides to cool down for 20mins. (Do not add cold tap water into the pressure cooker).


5. Monitor Workflow tab




- 5.1 While slides are being cooled down, make sure my slides tab under workflow tab will show the slides you entered.
 - 5.2 Determine and prepare required reagents using DAKO reagent bottles by registering them with Lot number and expiry date (5ml, 12ml and 25ml).
 - 5.2.1 Endogenous enzyme block (Dako, #SM801) or equivalent to 3% H2O2 (Sigma, #H1009)
 - 5.2.2 Primary antibody- MCM2(Cytosystem, #03/0912, 1in500)
 - 5.2.3 Labelled polymer (Dako, #SM802) or equivalent (DAKO anti-mouse HRP, #K4001)
 - 5.2.4 Substrate-Chromogen (SM803) or DAKO DAB substrate (Dako, #K3468) detection
 - 5.2.5 Counterstain-Harris Haematoxylin and Scotts tap water substitute
 - 5.2.6 Prepare wash buffer (DAKO, #DM831) or equivalent wash buffer (0.05% PBST) and deionised water and make sure waste bottle is empty
 - 5.3 Print and apply labels for each reagent's bottle and load it to reagent rack and the rack into the machine if necessary.
6. Wash slides in wash buffer for 5min
 7. Load slides into the black slide rack on the machine and moisten slides with wash buffer. (*Note:Make sure the slide rack is straight and even.)



8. Click the '**Instrument**' tab and choose 'AS1173D0903'
9. Prime the buffer and water. When the following message displays, click OK. In this step, please make sure in the buffer and water nozzles are not having any air bubbles. If there's then repeat priming the buffer and water.



10. Starting the run in the instruments tab
(Note: **Green** colour indicates normal operation message.
Yellow indicates instrument needs attention but will not spoil slides.
Conditions that cause yellow alerts include indications that reagents are required before the run can be processed.
Red indicates that you need to address the problem immediately.
Conditions that cause red alerts include system error messages.)




11. During the staining run, every 30mins keep checking the screen if there's yellow or red warnings. If so, fix the problems. Eg) Lack of reagents (yellow warning) – add more reagents)
12. End of the run, click the DONE button then the machine will purge the waste liquid.
13. Open the autostainer lid and remove completed slide racks and transfer each slide into designated wash buffer jar.
14. Dehydrate the slides in 50%, 80%, 100% and 100% alcohol for 30sec, 30sec, 2min and 2min respectively.
15. Clear the slides in xylene 3 x 5min
16. Mount slides with DPX.



5.0 Troubleshooting

PDF manual is available in Help



6.0 Personal protection -

A Howie laboratory coat and lab gloves must be worn at all times.

7.0 Spillages -

Always clean up any spills immediately after use, only you know what you have spilt and are aware of its hazard.

8.0 Training -

All staff should undergo training in this technique before performing procedure.

9.0 Related documents –

- 9.1 Risk assessments –COSHH RA 08 and 21
General RA 06 and 47
- 9.2 SOP SASoM/EQUIP/015
Use of the Pressure Cooker for Antigen Retrieval
SOP SASoM/METHOD/041
Haematoxylin and Eosin (H&E) Staining



9.0 Approval and sign off –

Author:

Name: In Hwa Um

Position: Post Doc

Signature: Date:

Management Approval:

Name: Peter Mullen

Position: Research Fellow

Signature: Date:

QA release by:

Name: Alex MacLellan

Position: QA Manager

Signature: Date:

