

Method Procedure

Document N	lumber:	SASoM/METHOD/141.v1
Title:	DNA & RNA	hydrolysis for analysis by mass spectrometry.
Version:	v1	
Author:	Jennifer Bre	)

Effective from:	24/08/2021		
Valid to:	24/08/2023		

SOP History		
Number	Date	Reason for Change
v1	24/08/21	Original

### 1.0 Purpose -

This SOP describes the current procedures for performing DNA & RNA hydrolysis for analysis by mass spectrometry in Laboratory 248/249 at the St Andrews School of Medicine (SASoM).

#### 2.0 Scope -

This SOP applies to all staff performing DNA & RNA hydrolysis for analysis by mass spectrometry in Laboratory 248 at the St Andrews School of Medicine (SASoM).

## 3.0 Responsibilities -

All staff performing DNA & RNA hydrolysis for analysis by mass spectrometry in this manner are responsible for ensuring that the methods are followed in accordance with this SOP. All staff must have read and signed the relevant risk assessment documents before performing this procedure.

Method Procedure



## 4.0 Procedure –

Reagents & buffers

Alkaline phosphatase #18011015 - Thermo DNAse I # 4716728001 - Roche Venom phosphodiesterase # P3134 - Sigma S1 nuclease # M5761 - Promega Calf intestinal phosphatase #M0525 - NEB LC-MS water: provided by CMC group (biology)

#### <u>DNA buffer</u>

 Tris HCl 20 mM pH7.9
 240mg

 NaCl 100 mM
 584 mg

 MgCl2 20 mM
 190.4 mg

 in 100 mL dH2O
 190.4 mg

#### <u>RNA buffer</u>

100 mM Tris HCl1.2 g50 mM NaCl292 mg10 mg MgCl295.2 mgIn 100 mL dH20Store both on bench

## For DNA hydrolysis

3.75 mL DNA buffer + 1  $\mu$ L alkaline phosphatase + 25  $\mu$ L DNAse I + 1.5  $\mu$ L venom Keep leftover in freezer for next batch of samples

#### For RNA hydrolysis

Per sample: 142 μL RNA buffer + 3.2 μL S1 nuclease + 5 μL CIP

#### <u>Protocol</u>

- 1. DNA and RNA to be collected according to the relevant SOPs as detailed in (i) SASoM-METHOD-129-DNA extraction using the QIAmp DNA Kit and (ii) SASoM-METHOD-093-Preparation of RNA using the Qiagen RNeasy kit.
- 2. Nanodrop for quantification

**N.B** Amounts required vary based on what target you are looking for (i.e 3  $\mu$ g DNA for dFdC but 10  $\mu$ g for FdUr) - talk with Alison and Jen.

- 3. Add the nucleic acid in a total volume of 10  $\mu$ L in LC-MS water (it is possible to go up to 50  $\mu$ L if the starting concentration is low to start with)
- 4. Add 150  $\mu$ L hydrolysis buffer
- 5. Briefly vortex
- 6. Incubate at 37°C overnight (16h min) in a thermocycler
- 7. Prepare samples spreadsheet for mass spec analysis



school of MEDICINE

Method Procedure

## 5.0 Personal protection -

A Howie coat must be worn at all times.

#### 6.0 Spillages –

Always clean up any spills immediately after use, only you know what you have spilt and are aware of its hazard. Spillages should be mopped up with paper towel, disinfected with 70% ethanol and finally washed with detergent.

#### 7.0 Training –

All staff should undergo training in this technique before performing the procedure.

#### 8.0 Related documents –

8.1 Risk assessments – RA12771: DNA & RNA hydrolysis for analysis by mass spectrometry

## 9.0 Approval and sign off –

Author:		
Name:	Jennifer Bre	
Position:	Post Doctorate	
Signature:	Jere -	Date: 06/09/21
Management Appr	oval:	
Name:	Peter Mullen	
Position:	SOP Administrator	
Signature:	Peter Muller	Date: 06/09/21
QA release by:		
Name:	John O'Connor	
Position:	QA Manager	
Signature:	Jr 5	Date: 07/09/2021



Method Procedure

# STANDARD OPERATING PROCEDURE

Please sign below to indicate you have read this S.O.P and understand the procedures involved.

NAME	POSITION HELD	SIGNATURE	DATE